

## Get started



### 1 Download the Neosensory app



- Experience Music mode
- Rest easy with Sleeping mode
- Adjust vibration strength
- Feel sample sounds




### 2 Put on your Buzz




Make sure it's snug. Then slide the pin to secure the end of the band.

### 3 Press and hold the power button



When your Buzz powers on, you will see three green lights.

### 4 Make any noise to feel your Buzz



### 5 Adjust to your environment



Switch between low, medium, and high microphone sensitivities.

### Additional info



Charge your Buzz in 30 min for a full day of use.



Buzz is water resistant but not waterproof. Wear in light rain but not in the shower.



For more instructions, download the printout at [www.neosensory.com/buzzinstructions](http://www.neosensory.com/buzzinstructions)

### Need more help?



Check out the Help Center in the Neosensory app or contact us at [www.neosensory.com/support](http://www.neosensory.com/support)